

Introducing a Breakthrough in PT/INR Testing

The Coag-Sense® Prothrombin Time (PT)/INR Professional Monitoring System

The First Portable System with Direct Clot Detection Technology

[Click Here](https://wilburnmedicalusa.com/coag-sense-ptinr) to Purchase or Visit <https://wilburnmedicalusa.com/coag-sense-ptinr>

Limited Time Promotion

Receive a FREE
Coag-Sense PT/INR System
with the Purchase of
3 Boxes of 50 Test Strips

\$570.00



Convenient,
in-office use



Simple sample
transfer tubes
included



Direct clot
detection
technology



Patients benefit
with more-
frequent testing



Control extended
from office to home

- First system with Direct Clot Detection Technology
- Accurate results in less than one minute
- Ability to run whole blood or plasma samples
- Not dependent on hemoglobin and hematocrit levels
- Superb precision (CV 2.5%)
- Small sample size (10 µl minimum) for easy dosing
- True functional controls using actual thromboplastin and plasma
- Barcoding on each test strip - no calibration code input or chip required



Visit or Call Us Today!

1-877-WILBURN (1-877-945-2876)

**146 Furlong Industrial Drive
Kernersville, NC 27284**

WilburnMedicalUSA.com



We also have a Coag-Sense Home User Promotion that you can Print or Download and E-mail for Your Patients. [Click Here](#)

or go to https://wilburnmedicalusa.com/content/pdf/Coag-Sense_Self-Test_Home_User_Promotion.pdf



The Coag-Sense PT/INR System:

ACCURATE PT/INR TESTING IS AT HAND WITH DIRECT CLOT DETECTION TECHNOLOGY

Our introductory offer includes:

- 3 boxes of test strips (quantity 150) plus 6 high and 6 low controls
- 1 analyzer
- 1 professional user manual
- 1 professional quick reference guide
- 1 bag of auto single-usage lancets (quantity 12)
- 4 AA alkaline batteries
- 1 A/C power adapter



Ordering Information

06F23-03 COAG-SENSE PROFESSIONAL PT/INR SYSTEM 3-BOX PROMO KIT

03P60-01 COAG-SENSE PROFESSIONAL PT/INR SYSTEM

03P56-50 COAG-SENSE PROFESSIONAL TEST STRIP KIT- 50

Click Here to Purchase or Visit <https://wilburnmedicalusa.com/coag-sense-ptinr>

© CoaguSense, Inc.
48377 Fremont Blvd., STE 113
Fremont, CA 94538
(866) 903-0890
(510) 226-6540 (Fax)
www.coagusage.com

Coag-Sense is a trademark of CoaguSense, Inc.

Cleared for professional use	Yes, CLIA waived
Cleared for home use	Yes
Sample size	1 drop (10 µL) capillary blood minimum
Sampling method	Single-use autolancet; micropipette transfer
Device specifications	Portable
Size	3 in (7.6 cm) x 6.5 in (16.5 cm) x 5.75 in (14.5 cm)
Weight	1.2 lb (0.5 kg) with batteries
Test strips packaging	50 Individually pouched, plus 2 low and 2 high controls (Avoids Contamination Issues)
Storage	35F°-72 F°, up to 24-month outdate; refrigeration not required
Coding	Barcoded, no code chip or entry required
Heparin sensitive	Yes
Clot detection principle	Direct; measures actual time required for clotting
INR range	0.8–8.0
Precision Ave. CV %	2.5% (Professional Use)
ISI	~1.0
Prothrombin time (PT)	8–80 seconds
Analysis time	<1 minute (Actual Prothrombin Time)
Hematocrit Range	15-60% (No RBCs Required)
QC	True plasma controls included with test strips
Memory	100 tests
Barcoded test strips	Yes, room temp. storage with 24 month outdating
Power source	4 AA batteries or AC adapter (professional system)
Warranty	1 year with immediate replacement
Maintenance	None required

For in vitro diagnostic use

To perform a CLIA-waived test, a certificate of waiver (COW) is required. For more information please contact your local distribution representative or CoaguSense, Inc. at **(866) 903-0890**

Professional Intended Use: The CoaguSense Prothrombin Time (PT)/INR Monitoring System is a diagnostic device that provides quantitative prothrombin time (PT) results, expressed in seconds and INR units. It uses fresh capillary whole blood. It is intended for use by health care professionals at the point of care to monitor patients who are on warfarin-type (coumarin) anticoagulation therapy. The device is not intended to be used for screening purposes.