



**YOUR TRUSTED  
LOW COST EQUIPMENT & SUPPLY  
RESOURCE SINCE 1997  
1-877-WILBURN (945-2876)**

## **Wilburn Medical is offering the following Infrared Digital Thermometer**

### **Product Code: IDT**

**Product Description:** 32-42.9°C. (89.6-109.2°F) (Accuracy below)  
**Response Time:** 1 second  
**Measure Distance:** 0-2"  
**Auto Shut Off Time:** 30 seconds  
**Small Sizes:** 5.9" x 1.5" x 1.7"  
**Certification:** CFDA, CE-mark, FDA 510(k) Exempt  
**Packaging:** 22.6" x 13.8" x 6.1",  
Net Weight: 5.9 lbs,  
Gross Weight: 13.7 lbs

**Measurement Range:** Human Body Temperature 32.0°C - 43.0°C  
**Display Resolution:** 0.1°C  
**Measurement Place:** Forehead  
**Accuracy:** 32.0° - 34.9°C, maximum permissible error =/- 0.3°C  
35.0°C - 42.0°C, maximum permissible error +/- 0.2°C  
42.1°C - 42.9°C, maximum permissible error +/- 0.3°C  
**Display Screen:** LCD screen with 4-digit and special icons  
**Operating Environment:** 10°C - 40°C  
**Max. Relative Humidity:** <85%  
**Storage Environment:** -20°C - +55°C  
**Battery:** 2 x 1.5V AAA batteries



**CALL 1-877-WILBURN (945-2876)  
OR VISIT WILBURNMEDICALUSA.COM TO ORDER!  
140 Furlong Industrial Dr • Kernersville, NC 27284**

