



JXB-183

Non-Contact Infrared Thermometer

HELPS PREVENT
THE SPREAD OF
INFECTIOUS DISEASE

Provides quick & accurate measurement
of core body temperature

- No probe covers required - significant annual savings
- No contact with resident's skin reduces the risk of cross-contamination
- Ideal for sleeping or combative residents
- FDA cleared, clinical-grade thermometer for multi-patient and individual use
- Ideal for F-Tag 257 compliance (comfortable & safe temperature levels)

FEATURES

- Non-contact reading in one second from 1.2 - 2 inches distance from forehead
- LCD digital display screen with high temperature alarm
- Stores the previous 32 temperature measurements
- Temperature can be displayed in either °F or °C
- Automatic power-off (30 seconds) to conserve energy
- Modes for body, surface and room (ambient) temperature
- Visual/audible alarm if temperature exceeds 100.4°F
- Long battery life
- 1 Year limited warranty - manufacture defect only

ACCURATE
and precise

INSTANT
1 second
measurement

AUDIBLE
option



**CALL 1-877-WILBURN (945-2876) OR
VISIT WILBURNMEDICALUSA.COM TO ORDER!**

VALIDATED ACCURACY

The JXB-183 non-contact thermometer has a validated accuracy of +/- 0.2 °C (0.4 °F) over the clinical range of valid temperature.*

* Standard for ISO 13485