



JXB-183

Non-Contact Infrared Thermometer

Device Type	Non-contact infrared thermometer
Measurement Distance	1.2" to 2" (3 to 5 cm) from forehead
Measuring Range	Body temperature: 89.6°F to 109.2°F (32°C to 42.9°C) Surface temperature: 32°F to 140°F (0°C to 60°C) Ambient temperature: 50°F to 104°F (10°C to 40°C)
Displayed Temperature Range	89.6°F to 109.2°F (32°C to 42.9°C)
Display Type	Large, 2-color LCD with highly accurate 1/100 degree reading
Display Resolution	0.1°F (0.1°C)
Accuracy	89.6°F to 96.62°F=±0.6°F (32°C to 35.9°C=±0.3°C) 96.8°F to 102.2°F=±0.4°F (36°C to 39°C=±0.2°C) 102.3°F to 109.2°F=±0.6°F (39.1°C to 42.9°C=±0.3°C)
Operating Temperature	50°F to 104°F (10°C to 40°C)
Operating Humidity	Up to 85% non condensing
Storage Temperature	-4°F to 131°F (-20°C to 55°C)
Storage Humidity	Up to 95% non condensing
Consumption	≤300 mW
Calibration	Preset at factory. Device re-calibrates each time device is turned ON.
Memory	Stores up to 32 readings
Automatic Power Off	30 seconds after last measurement
Batteries	2 AAA batteries
Unit Weight	75g without batteries
Warranty	1 Year limited warranty



**CALL 1-877-WILBURN (945-2876) OR
VISIT WILBURNMEDICALUSA.COM TO ORDER!**