

The CardioChek® PA Analyzer

The CardioChek PA is a cost-efficient, user-friendly diagnostic point-of-care analyzer, which accurately determines multiple patient values.

Features & Benefits



FAST

Results in less than 2 minutes



PORTABLE

Handheld, battery-operated



PATIENT-FRIENDLY

Small sample (fingerstick) 15-40 µL



QUALITY

CLIA-waived, FDA-cleared, CE-marked;
CRMLN-certified for CHOL and HDL test strips;
Meets National Cholesterol Education
Program (NCEP) Guidelines



HIGH RETURN ON INVESTMENT (ROI)

Rapid throughput and low waste,
with minimal start-up cost



Available PTS Panels® Test Strips

Customize your diagnostic offerings with a variety of PTS Panels test strips for use with the CardioChek® PA analyzer. Our test strips offer an easy, efficient method for point-of-care diagnostic testing and Health and Wellness Screening.



Multi-Analyte



LIPID Panel Test Strip

Measures: Total Cholesterol, HDL Cholesterol & Triglycerides
Calculates: LDL Cholesterol, Total Cholesterol/HDL Ratio



CHOL+HDL+GLU Panel Test Strip

Measures: Total Cholesterol, HDL Cholesterol & Glucose
Calculates: Total Cholesterol/HDL Cholesterol Ratio



METABOLIC CHEMISTRY Panel Test Strip

Measures: HDL Cholesterol, Triglycerides & Glucose

Dual Panel



CHOL+GLU Panel Test Strip

Measures: Total Cholesterol & Glucose



CHOL+HDL Panel Test Strip

Measures: Total Cholesterol & HDL Cholesterol
Calculates: Total Cholesterol/HDL Cholesterol Ratio

Single Test Strips



TOTAL CHOLESTEROL Test Strip

Measures: Total Cholesterol



HDL CHOLESTEROL Test Strip

Measures: HDL Cholesterol



TRIGLYCERIDES Test Strip

Measures: Triglycerides



GLUCOSE Test Strip

Mesure : glucose

Other test strips available include PTS Panels LDL Cholesterol (direct, not calculated), Creatinine (special order only), and Ketone (special order only).



PTS Diagnostics

7736 Zionsville Road, Indianapolis, IN 46268 USA
Office: +1 (317) 870.5610 FAX: +1 (317) 870.5608
Toll-Free in the U.S. : +1 (877) 870.5608

ptsdiagnostics.com