



JMS Co., Ltd.

12-17 KAKO-MACHI, NAKA-KU  
HIROSHIMA 730-8652 JAPAN  
PHONE: +81-(0)82-243-5505  
FAX: +81-(0)82-243-1976

No. 230204-C

Date: February 24, 2023

## Certificate

JMS DISPOSABLE SYRINGE 551721050,

This is to certify that the 100 mL syringe family is dose audit tested via Method 1 of ISO 11137-2 on a quarterly. All tests are found to be acceptable to prove a  $10^{-6}$  sterility assurance for the listed products when processed via EB radiation at a dose of 20-55kGy.

JMS Co., Ltd.

Makoto Goto

Head of Quality Assurance Section