

WUSA/100

Single Use Sterile Luer Lock Syringe

General Information		Materials	
Volume	100 mL (rated) marked to 120 mL	Barrel	Polypropylene (PP)
Scale Marks	Marked in 1 mL increments	Plunger	Polypropylene (PP) Polycarbonate (PC)
Manufacturer & REF Code	JMS Ltd., Hiroshima Japan Item# 551721050 ISO 13485:2016 Certified Certificate# Q5 024761 0166 Rev. 05	Gasket	Thermo Plastic Elastomer (TPE)
		Lubricant	Silicone
		Package (Tear Strip)	Paper
Classification	Class II / FDA 510K# K132321	Package (Pouch)	Polyethylene (LDPE)
Dead Space	Less than 0.2 mL		
Sterilization	EB Radiation at 20- 55kGy		
Shelf Life	3 years from DOM		
Melamine Free Latex Free BPA Free Phthalate Free Biocompatibility - Tested to ISO 10993 Extractables - Tested to ISO 7886-1			