



**Wilburn Medical USA**  
**1-877-WILBURN (1-877-945-2876)**

Wilburn Medical Inc.  
140 Furlong Industrial Drive  
Kernersville, NC 27284  
Phone: 1-336-996-6820 Fax: 1-336-992-0847  
[wusainfo@wilburnmedical.com](mailto:wusainfo@wilburnmedical.com)

[WilburnMedicalUSA.com](http://WilburnMedicalUSA.com)

**PERSONAL PROTECTION**

# Complete Face Protection

## Opti-Shield

Offers full face protection from splashes, sprays, splatters, or droplets of potentially infectious materials that pose a hazard to the eyes, nose, or mouth. The wraparound design of the Opti-Shield allows complete, distortion free viewing. The foam band keeps the shield a comfortable distance from the face and does not interfere with most eyeglasses.

### Product Features

- Adjustable foam headband
- Lies flat against chest when not in use
- Anti-fog shield
- Lightweight
- Cost effective

### Your Benefits

- Guarantees a perfect fit for everyone
- Stays out of the way
- Gives optimum visibility while wearing
- All-day comfort
- Will last as long as more expensive visors

### Comfort by Design

Opti-Visor gives complete facial protection and offers unsurpassed comfort. It is easy to both don and remove and fits over glasses. Our unbreakable, fog-free shields are easily adjustable and can be worn around the neck when not in use. At less than one ounce, Opti-Shield is comfortable even for long periods of wear.



Item Number: WUSA/SOURCE-70952  
50 Shields per Case

